

## **Protocol: Extraction of human A $\beta$ from APP Tg mouse brains into TBS-soluble, guanidine HCl-soluble fractions**

### **TBS-soluble fraction**

1) Make 10 ml of fresh TBS extraction buffer (sterile 1x TBS + 100  $\mu$ l of 200 mM 1,10-phenanthroline + 100  $\mu$ l of 500 mM EDTA + 1 tab of Roche complete EDTA-free Mini protease inhibitor).

#### TBS extraction buffer final concentrations:

- 140 mM NaCl
- 3 mM KCl
- 25 mM Tris-HCl, pH 7.4
- 1 tab/10 ml of Roche complete EDTA-free Mini protease inhibitor
- 5 mM EDTA
- 2 mM 1,10-phenanthroline

- 2) Weigh frozen hemisphere quickly, before it thaws. After weighing, keep brain in a weigh boat on dry ice until homogenization.
- 3) Put 4 brain volumes (4  $\mu$ l/mg of brain) of ice-cold TBS extraction buffer into a 2 ml Potter-Elvehjem (PE) homogenizer on ice, then add the brain hemisphere.
- 4) Attach the pestle to a motorized rotor and homogenize sample. Do 20 strokes.
- 5) Transfer homogenate (with a long-necked transfer pipette) to an ultracentrifuge tube, and store on ice until all samples are homogenized.
- 6) Weigh tubes and arrange them into weight-matched pairs. Put each member of the pair across from its match in the rotor.
- 7) Spin at 100,000xg for 1 hr at 4°C. Save the supernatant as the “TBS-soluble fraction”.

### **Guanidine HCl-soluble fraction**

- 1) Add 4 volumes (based on initial hemisphere weight) of cold 6.25 M guanidine HCl in 50 mM Tris, pH 8.0 to the pellet.
- 2) Vortex the tube until the pellet becomes dislodged from the bottom.
- 4) Transfer sample to a clean PE homogenizer. Make sure all pellet fragments are transferred and not stuck in the tube or pipette.
- 5) Using the motorized pestle to homogenize sample with 20 strokes.
- 6) Transfer the homogenate with a long-necked transfer pipette to a new tube.

7) Nutate for ~2 hr at RT, then spin at 4°C 20,800xg for 20 minutes. Save supernate as the “GuHCl-soluble fraction,” and discard the small viscous pellet.