

## Neprilysin Activity Assay

Substrate DAGNPG (N-dansyl-Ala-Gly-D-nitro-Phe-Gly, modeled after enkephalins with an aromatic moiety in the P'1 position and a short residue in the P'2 position) is an internally quenched substrate of NEP. NEP degrades DAGNPG to DAG ( $K_m = 45 \mu\text{M}$ ,  $V = 0.65 \mu\text{mol/mg protein/min}$ ), and the liberated dansyl group can be excited at a wavelength of 342 nm and emits at 562 nm. ACE can also degrade DAGNPG, therefore an ACE inhibitor is necessary to resolve NEP activity (must pre-incubate for 10 min, using at least  $0.5 \mu\text{M}$  captopril). NEP has optimal activity in 50 mM Tris HCl, pH 7.4, and inhibited by 20 nM Thiorphan ( $K_i = 3 \text{ nM}$ ).

Prepare a stock 25 mM DAGNPG solution in methanol (store at  $4^\circ\text{C}$ ).

1. Prepare 100 $\mu\text{g}$  lysate in 50 mM Tris HCl, pH 7.4 with captopril. The final volume (per reaction, including substrate) should be 200  $\mu\text{L}$ .
2. Make a 1 mM DAGNPG solution in Tris. Add DAGNPG to the lysate for a 50  $\mu\text{M}$  final concentration, mix and incubate at  $37^\circ\text{C}$  for desired time points.
3. Stop the reaction by heating to  $100^\circ\text{C}$  for 5 min. Store tubes on ice until all time points are reached.
4. Spin 5,000g x 5 min at RT.
5. Remove 175  $\mu\text{L}$  sample and dilute into 400  $\mu\text{L}$  Tris. Vortex to mix and add 200  $\mu\text{L}$  to a 96-well plate.
4. Read on the Victor2 multilabel plate reader (excitation, 355; emission, 550, optimal wavelengths are 342/562).